

## ASSAY FOR THE SECRETED ALKALINE PHOSPHATASE (SEAP)

- 1> - From the 48-hour cotransfected cell cultures, remove 250  $\mu$ l of each culture supernatant and transfer into Eppendorf tubes. To be on the safe side, maintain the cultures at 37°C until satisfactory data have been obtained.
- 2> - Heat samples at 65°C for 5 minutes to inactivate endogenous phosphatases (Seap is stable at 65°C). *homoarginine also inactivate endogenous phosphates*
- 3> - Centrifuge for 2 minutes at room temperature in a Microfuge.
- 4> - Transfer <sup>100  $\mu$ l</sup> supernatants to new Eppendorf tubes. At this stage, samples may be stored at -20°C indefinitely.
- 5> - In an Eppendorf tube, add 100  $\mu$ l of 2 x Seap buffer to 100  $\mu$ l of sample. As a Zero standard make up a mixture in triplicate substituting sample with water. Mix on a Vortex.
- 6> - Transfer the contents of each tube to a well of a flat bottom microtiterplate. Avoid creating air bubbles.
- 7> - Incubate plate at 37°C for 10 - 15 minutes.
- 8> - During this incubation make up the p-nitrophenylphosphate solution (Seap enzyme substrate) and prewarm it at 37°C for 5 minutes.
- 9>- Add 20  $\mu$ l of the substrate solution to each well, preferably using a multipipetter.
- 10>- Using an ELISA microplate reader with an automatic shaker and incubator unit, measure the OD at 405 nm at regular intervals (e.g. every 5 minutes) over 60 minutes while the plate is being incubated at 37°C. (Program a 5-second shaking before any reading).
- 11>- Calculate the levels of Seap activity at a point on the curve when the changes in OD are linear with respect to time (e.g. at 15 - 30 minutes).

\*\*\*\*\* \* \*\*\*\*\*

## Buffer and Chemicals

## ✓ 2 x Seap Buffer (for 50 ml)

Amount	Stock	Final Conc.
10.51 g*	diethanolamine (100% sol.)	2 M
50 $\mu$ l	1 M MgCl <sub>2</sub>	1 mM
226 mg	L-homoarginine	20 mM

(\*) Weigh exactly 10.51 g of diethanolamine in a tared beacker. Add distilled water up to 45 ml. Stir to homogenize. Add 50  $\mu$ l of 1 M MgCl<sub>2</sub> while stirring. In a separate 15-ml conical tube, dissolve 226 mg of L-homoarginine in 2-3 ml of distilled water. Add the solution to diethanolamine/MgCl<sub>2</sub> solution under constant stirring. Bring up to 50 ml.

120 mM p-nitrophenylphosphate is made in 1 x Seap buffer (fresh).

Amount G (mg) = (120 mM x 371.12 x vol)/1000

Where: vol = [(# wells x 20  $\mu$ l) + 100  $\mu$ l extra]/1000

\* Make the solution in 1 x Seap buffer (make fresh)

For 51 wells of the microplate (48 samples + 3 blanks), dissolve 50 mg p-nitrophenylphosphate in 1.120 ml of 1 x Seap buffer.

→ 600  $\mu$ l stock + 600  $\mu$ l H<sub>2</sub>O

## Chemicals

Name	Cat No	Storage
Diethanolamine (Fluka)	31589	Room T
L-homoarginine hydrochloride (Sigma)	H-1007	4°C
p-nitrophenylphosphate* (Fluka)	71768	4°C

(\*) Also known as 4-nitrophenylphosphate disodium salt hexahydrate

0.21806

0.226

# SEAP ASSAY SHEET

## I - Assay Title:

\* Assay #: 72

Date #: [REDACTED]

Investigator's Name: Wen

\*Test Compounds: Mal 4 & Mal 43Na<sup>+</sup>

transfected at [REDACTED] 11:30

\*Concentrations:

## II - DNA Transfection: Ratio (2:1) HIV/SEAP:pcTAT

\*HIV/Seap: 0.4 µg/well x 40 = 16.0 µg ==> From 0.376% stock: 42.55 µl

\*pcTAT: 0.2 µg/well x 40 = 8.0 µg ==> From 0.639% stock: 12.52 µl

\*Total DNA (µg) = 24 µg

\* Vol. Cellfectin = Total DNA x 6 = 144 µl

\* Vol. 150 mM NaCl = Cellfectin (µl) / 0.6 = 240 µl

HIV/SEAP  
42.55 µl

pcTAT  
12.52 µl

(A) 240 µl SALT

(B) Cellfectin 144 µl

\* Transfection cocktail (µl)/ well:

$$\frac{240 + 240 + 144 + 42.55 + 12.52}{40} = 16.98 \mu l$$

## III. Linbro® 24 flat bottom well of 17 mm

A

10	10	10	10	10	10
20	20	20	20	20	20

B

10	10	10	10	10	10
20	20	20	20	20	20

## Preparation of Drug conc.

(NOTE)

C

C	C	C	C	C	C
CT	CT	CT	CT	CT	CT
C DMSO	C DMSO	C DMSO	C DMSO	C DMSO	C DMSO
CT DMSO	CT DMSO	CT DMSO	CT DMSO	CT DMSO	CT DMSO

40: 78  
60: 97

40: 23  
60: 46

REPEATED-READS #: 0025

07:11:06.56 PM

TABLE OF ABSORBANCE VALUES

FILE: FILE 3 TITLE: FILE 3 REPORT

1 2 3 4 5 6 7 8 9 10 11 12

Mal. 4 DMSO											
A	0#000-0#002	0#004-0.217	-0.215	-0.216	0.064	0.080	0.089	0.126	0.133	0.146	
					0.078			0.135			
B	-0.210	-0.214	-0.205	-0.215	-0.213	0.151	0.153	0.207	0.124	0.153	0.155
						0.170			0.144		
C	0.070	0.016	0.060	0.044	-0.214	-0.214	0.105	0.089	0.131	0.171	0.169
			0.048				0.108			0.170	
D	0.689	0.649	0.663	0.756	-0.208	-0.210	0.126	0.122	0.120	0.172	0.197
			0.689				0.123			0.185	
C DMSO											
E	0.094	0.100	0.081	0.110	-0.212	-0.209	0.093	0.067	0.087	0.101	0.100
			0.096				0.082			0.094	
F	0.680	0.711	0.680	0.614	-0.211	-0.213	0.830	0.776	0.832	0.649	0.587
			0.671				0.813			0.579	
C NaOH											
G	0.064	0.067	0.080	0.079	-0.203	-0.211	0.089	0.114	0.100	0.114	0.118
			0.073				0.101			0.110	
H	0.784	0.708	0.643	0.684	-0.209	-0.205	0.681	0.720	0.609	0.481	0.437
			0.705				0.670			0.452	

Mal. 4 - DMSO Ref: DMSO  
% inhibition

10  $\mu\text{g/mL}$ : 84  
20  $\mu\text{g/mL}$ : 97  
40  $\mu\text{g/mL}$ : 98  
60  $\mu\text{g/mL}$ : 97

Mal. 4 - NaOH Ref: NaOH  
% inhibition

10  $\mu\text{g/mL}$ : 0 (-16)  
20  $\mu\text{g/mL}$ : 10  
40  $\mu\text{g/mL}$ : 23  
60  $\mu\text{g/mL}$ : 46

## SEAP assay of NDGA derivatives

1. Detach COS7 cells ~~from~~ <sup>adding</sup> four 90mm culture plates by 0.05% Trypsin (0.05% Trypsin in CMF-PBS, 1mM EGTA) 0.5ml after 2 times CMF-PBS Wash.
2. Inactivate trypsin by adding 5ml complete <sup>IMDM</sup> medium (10% Fetal Bovine Serum, antibiotics) to each plate.
3. Suspend the cells by gentle pipetting. Count by Hemocytometer  $160 \times 10^4$  cells/cm<sup>3</sup> is the cell density.
4. Seed ~~the~~ 80ul cell suspension into 10 Linbro 24-wells culture plates which <sup>were</sup> precoated with 0.1% gelatin & contained 0.3ml complete medium each well. ( $1.3 \times 10^5$  cells/well) (18:00)
5. Transfect the cells after 26h incubation (37°C, 95% Air-5% CO<sub>2</sub>) by Adding 30ul ppt (0.6-ug DNA/well) to each well that containing 0.3ml fresh complete medium. The control wells accept same ppt without DNA. (20:00)  
Ca-PO<sub>4</sub> ppt preparation:
  - a. for control wells: Ca<sup>2+</sup> soln (0.25M CaCl<sub>2</sub>, 0.025M HEPES, pH 7.1) 1ml drop into 1ml bubbling PO<sub>4</sub><sup>2-</sup> soln (0.28M NaCl, 0.025M HEPES, ~~0.0015~~ <sup>0.0015</sup> M Na<sub>2</sub>HPO<sub>4</sub>, pH 7.1).
  - b. for transfect wells:  
Three 2-ml tubes containing ~~pBC12/HIV/SEAP 18ug, pBC12/CMV/tz 9ug~~ 675ul PO<sub>4</sub><sup>2-</sup> soln were added Ca<sup>2+</sup> soln 675ul ~~con~~ with DNA, (pBC12/HIV/SEAP 18ug, pBC12/CMV/tz 9ug) dropwise with bubbling. then sit 30min before use.
6. Incubate at 37°C for 18h.

7  
Map of Plates:

I

	C			C #1 20μM
	C Mal.4 20μM			C #2 20μM
	CT			CT #1 20μM
	CT Mal.4 20μM			CT #2 20μM

A # Mal.4

C <del>0μM</del>	C 10μM	C 60μM
C Mal.4 3μM	C 30μM	C 100μM
CT 0μM	CT 10μM	CT 60μM
CT 3μM	CT 30μM	CT 100μM

C #2

C 0μM	C 10μM	C 60μM
C 3μM	C 30μM	C 100μM
CT 0μM	CT 10μM	CT 60μM
CT 3μM	CT 30μM	CT 100μM

E #4

C 0μM	C 10μM	C 60μM
C 3μM	C 30μM	C 100μM
CT 0μM	CT 10μM	CT 60μM
CT 3μM	CT 30μM	CT 100μM

G #6

C 0μM	C 10μM	C 60μM
C 30μM	C 30μM	C 100μM
CT 0μM	CT 10μM	CT 60μM
CT 30μM	CT 30μM	CT 100μM

II

	C #3 20μM	C #4 20μM	C #5 20μM	C #6 20μM
	C #7 20μM	CT #4 20μM	CT #5 20μM	CT #6 20μM
	CT #3 20μM			
	CT #7 20μM			

B #1

C 0μM	C 10μM	C 60μM
C 3μM	C 30μM	C 100μM
CT 0μM	CT 10μM	CT 60μM
CT 3μM	CT 30μM	CT 100μM

D #3

C 0μM	C 10μM	C 60μM
C 3μM	C 30μM	C 100μM
CT 0μM	CT 10μM	CT 60μM
CT 3μM	CT 30μM	CT 100μM

F #5

C 0μM	C 10μM	C 60μM
C 3μM	C 30μM	C 100μM
CT 0μM	CT 10μM	CT 60μM
CT 3μM	CT 30μM	CT 100μM

H #7

C 0μM	C 10μM	C 60μM
C 3μM	C 30μM	C 100μM
CT 0μM	CT 10μM	CT 60μM
CT 3μM	CT 30μM	CT 100μM

C = Control wells (without DNA)  
CT = Transfect wells (with DNA)

8. Remove growth medium, add 0.500 ml medium containing test compounds. The final concentration of DMSO is 0.3% ( [REDACTED] )

Δ Dilution of test compounds:

Stock: Mono Me: 10.53  $\mu\text{g}/\mu\text{l}$  Use 3  $\mu\text{l}$  in 1 ml medium = 100  $\mu\text{M}$

DiMe: 11.0  $\mu\text{g}/\mu\text{l}$  "

TriMe: 11.41  $\mu\text{g}/\mu\text{l}$  "

tetraMe: 11.93  $\mu\text{g}/\mu\text{l}$  "

For lower concentration, dilute the stock by DMSO, & use 3  $\mu\text{l}$  in 1 ml medium to keep the final DMSO concentration is 0.3%.

Δ DMSO stock was ~~added~~<sup>mixed</sup> into medium just prior to use by vortexing.

9. Incubate for 48h. Remove 200  $\mu\text{l}$  medium from each well. for SEAP assay. ( [REDACTED] 14:00 ).

10. The medium used in SEAP assay is 10  $\mu\text{l}$  for each sample.



REPEATED READS #1 0025

05:50:03.28 PM

# TABLE OF ABSORBANCE VALUES

FILE: FILE 3 TITLE: FILE 3 REPORT

1 2 3 4 5 6 7 8 9 10 11 12

I							II						
A	-0#024	0#005	0#020	-0.217	-0.216	-0.217	0.011	0.012	0.013	0.004	0.005	0.002	
							0.012			0.004	0.007	0.001	
E	0.012	-0.008	0.012	0.012	0.020	0.008	0.002	0.004	0.005	0.003	0.019	0.003	
							0.004						
C	-0.004	0.004	0.016	0.004	0.008	-0.001	0.585	0.704	0.631	0.558	0.394	0.421	
							0.146	#3		0.530	0.411	0.416	
D	1.126	1.040	1.102	0.801	0.766	0.811	0.347	0.373	0.326	0.502	0.428	0.410	
Mark							0.349	#4		#4	#5	#6	
E	1.010	0.855	0.915	0.890	0.753	0.799	0.018	0.020	0.007	0.012	-0.006	0.023	
							0.019			0.010		0.009	
F	-0.213	-0.216	-0.215	-0.213	-0.213	-0.216	0.008	0.014	0.007	-0.010	0.010	0.036	
							0.011		-0.002			0.023	
G	0.007	0.011	0.013	0.005	0.001	0.017	1.323	1.433	1.105	1.111	0.292	0.269	
							1.378		1.108		0.281		
H	0.007	-0.001	0.016	0.003	0.009	0.026	1.284	1.234	0.783	0.715	0.165	0.162	
							1.559		0.749		0.164		

Mal. % inhibition

3 μM : 8.2

10 μM : 19.2

30 μM : 44.7

60 μM : 80.0

100 μM : 89.6

07:46:07.64 PM

TABLE OF ABSORBANCE VALUES

FILE: FILE 3 TITLE: FILE 3 REPORT

1 2 3 4 5 6 7 8 9 10 11 12

A	0#002-0#004 0#003-0.211-0.209-0.212 0.008 0.016 0.003-0.001 0.002 0.016										
	0.012 0.001 0.009										
B	0.002 0.002 0.012 0.012 0.015 0.008 0.003 0.009 0.000-0.005 0.010 0.009										
	0.002 0.012 0.011 0.006 -0.003 0.010										
C	-0.005 0.003 0.008 0.004 0.007-0.001 1.200 1.149 0.871 1.088 0.374 0.357										
	0.001 0.004 0.003 1.175 0.980 0.366										
D	1.385 1.306 1.025 1.098 0.423 0.464 1.067 1.185 0.570 0.530 0.236 0.228										
	1.346 1.062 0.444 1.126 0.550 0.231										
E	1.288 1.327 0.905 1.005 0.330 0.314 0.017 0.012 0.002 0.007 0.003 0.016										
	1.308 0.955 0.322 0.015 -0.005 0.011										
F	-0.202-0.210-0.210-0.210-0.209-0.209 0.006 0.010 0.009 0.002 0.004 0.010										
	0.008 0.006 0.010										
G	1.239 1.354 1.209 1.239 0.420 0.350 1.122 1.205 0.921 0.933 0.208 0.218										
	1.291 1.224 0.385 1.164 0.927 0.213										
H	1.387 1.279 0.682 0.580 0.203 0.193 1.178 1.165 0.475 0.435 0.085 0.114										
	1.333 0.631 0.198 1.172 0.455 0.100										

#1	#2	#3	#4
% inhibition	% inhibition	% inhibition	% inhibition
3 $\mu$ M : -3.3(0)	3 $\mu$ M : 2.8	3 $\mu$ M : 3.7	3 $\mu$ M : -1.3(0)
10 " : 5.7	10 " : 21.9	10 $\mu$ M : 15.8	10 $\mu$ M : 19.8
30 " : 52.8	30 " : 29.4	30 $\mu$ M : 52.5	30 $\mu$ M : 60.9
60 " : 70.8	60 " : 67.8	60 $\mu$ M : 69.3	60 $\mu$ M : 82.4
100 " : 86.0	100 " : 76.3	100 $\mu$ M : 81.0	100 $\mu$ M : 92.2



REPEATED-READS #: 0025

09:06:41.21 PM

TABLE OF ABSORBANCE VALUES

FILE: FILE 3 TITLE: FILE 3 REPORT

1 2 3 4 5 6 7 8 9 10 11 12

A 0#000-0#008 0#010-0.210-0.208-0.210-0.209-0.209-0.210-0.208-0.209-0.209

B -0.207-0.208-0.209-0.208-0.208-0.209-0.208-0.209-0.209-0.208-0.209-0.211

G #6

H #7

C	0.012	0.014	0.012	0.022	0.017	0.020	0.011	0.014	0.011	0.014	0.019	0.015
	0.013		0.013		0.019		0.013		0.013		0.013	
D	0.007	0.013	0.015	0.014	0.026	0.012	0.010	0.019	0.004	0.020	0.038	0.043
	0.010		0.015		0.019		0.015		0.012		0.041	
E	1.192	1.134	0.699	0.687	0.277	0.262	1.391	1.489	0.658	0.623	0.317	0.300
	1.163		0.693		0.270		1.440		0.641		0.309	
F	1.215	1.138	0.374	0.400	0.202	0.184	1.163	0.974	0.380	0.317	0.189	0.239
	1.177		0.387		0.193		1.069		0.349		0.214	

G -0.205-0.208-0.208-0.206-0.199-0.208-0.206-0.207-0.207-0.207-0.209-0.210

H -0.207-0.207-0.204-0.206-0.206-0.204-0.189-0.207-0.207-0.206-0.206-0.210

#6

% inhibition

#7

% inhibition

~~200µM~~ 3µM : -1.5  
10µM : 41.0  
30µM : 67.7  
60µM : 78.2  
100µM : 84.9

3µM : 26.1  
10µM : 56.0  
30µM : 76.4  
60µM : 79.5  
100µM : 87.9

REPEATED-READS #: 0025

10:33:14.31 PM

TABLE OF ABSORBANCE VALUES

FILE: FILE 3 TITLE: FILE 3 REPORT

1 2 3 4 5 6 7 8 9 10 11 12

A -0#005-0#012 0#018-0.209-0.210-0.211-0.210-0.209-0.208-0.209-0.210-0.211

E -0.208-0.209-0.209-0.209-0.209-0.209-0.209-0.209-0.209-0.209-0.207-0.211

0.028 0.031 0.037 0.037 0.040 0.041 0.206-0.207-0.209-0.208-0.208-0.209  
0.030 0.037 0.041

D 0.015 0.028 0.036 0.028 0.047 0.035 0.208-0.210-0.208-0.208-0.208-0.209  
0.022 0.032 0.041

E 1.500 1.394 0.701 0.680 0.318 0.291 0.207-0.207-0.208-0.207-0.209-0.209  
1.447 0.691 0.305

F 1.131 1.033 0.385 0.417 0.222 0.185 0.207-0.208-0.207-0.206-0.208-0.212  
1.082 0.401 0.504

G -0.206-0.206-0.208-0.207-0.206-0.208-0.206-0.205-0.207-0.207-0.208-0.208

H -0.207-0.204-0.209-0.208-0.207-0.207-0.207-0.207-0.207-0.207-0.207-0.207

#5 % inhibition

3 μM : 25.2

10 μM : ~~53.8~~ 53.8

30 μM : 74.0

60 μM : 81.4

100 μM : 88.5

## 7 Inhibition of NDGA Derivatives on SEAP assay.

	Max. 4	#1	#2	#3	#4	#5	#6	#7
3 $\mu$ M	8.2	0.33	2.8	3.7	0.13	25.2	0.15	26.1
10 $\mu$ M	19.2	5.7	21.9	15.8	19.8	53.8	40.1	56.0
30 $\mu$ M	44.7	<del>52.8</del> 52.8	29.4	52.5	60.9	74.0	67.7	76.4
60 $\mu$ M	80.0	<del>52.8</del> 70.8	67.8	69.3	82.4	81.4	78.2	79.5
100 $\mu$ M	89.6	86.0	76.3	81.0	92.2	88.5	84.9	87.9

Figure B(1)




Figure B(2)

